**Congratulations you now own
a state of the art stethoscope
and thermometer made with
advanced technology and fine
materials.**

—

Contents

—

(English)

3	Indication for use, Getting Started
4-5	Thermometer setup
6-7	Taking Recordings - Thermometer
8	Stethoscope setup
9-10	Taking recordings - Stethoscope
11	Support, Standards
12	Troubleshooting
13	Warnings
14	Limited Warranty, Instructions
15	Limitations
16-17	Technical information
18	Symbols used

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Indication for use

Stethoscope (SPL999)

The CliniCloud stethoscope is a digital stethoscope intended for periodic recording of lung and heart sounds to a smartphone device, via the CliniCloud app. It is intended to be used for persons of all ages.

Thermometer (SPL 1024)

The CliniCloud non-contact forehead thermometer is an infrared thermometer intended for the periodic measurement of human body temperature for persons of all ages, to be used in conjunction with the CliniCloud app.

Getting Started

Download the CliniCloud iOS or Android app

<http://clinicloud.com/app>

App features

- Take recordings from the CliniCloud thermometer and stethoscope
- Save, browse, play, send and delete recordings
- Create user profiles for each member of your family
- Share your recordings with your doctor* or Doctor On Demand**

Notes:

Please check with your physician to see if they are amenable to receiving recordings via email.


US and Canadian customers only, for information on similar telemedicine services in your country please check

<http://clinicloud.com/partners>

What's included


- CliniCloud Non-Contact Thermometer
- CliniCloud Stethoscope
- Micro-USB to 3.5mm cord
- Quick Start Guide
- Coin cell battery (CR2032)
- 2 x AAA batteries

Thermometer setup

Thermometer setup

- a. Slide battery cover off
- b. Insert 2 x AAA batteries
- c. Slide Battery cover on

Taking Recordings – Thermometer


Step 1) Ensure Bluetooth function on your smartphone is turned on. Follow on-screen instructions to connect to the thermometer. The LED will turn on.

Step 2) Initiate recordings by following on-screen instructions in the app.


Step 3) Swipe thermometer across forehead from temple to temple, with the sensor facing the skin. Keep the device 0 - 5 cm (0 -2 inches) away from skin. The recording will take 3 seconds from initiation. The front LED will be flashing intermittently during recording.

Caution: If forehead is moist, please wipe dry before use. Ensure hair is out of the way of the sensor while taking the temperature. For optimum readings, allow at least 15 minutes to adjust to new room temperature conditions before use.


Note: Thermometer will automatically turn off after 2 minutes of inactivity


Step 1.

Step 2.

Step 3.




Information about the thermometer

The CliniCloud thermometer uses a medical grade infrared sensor to capture heat from the forehead. It takes multiple readings while scanning the forehead and reports the maximum temperature (usually measured over the temple area). The reading displayed is equivalent to oral thermometer temperature (and takes into account the environmental temperature of the room where the measurement is performed). A guide to interpreting the displayed temperature is provided below. Please note, normal temperature is a range which varies from person to person and fluctuates during the day.

Measurement	Range
35.8 - 37.4 °C (96.4 - 99.4 °F)	Normal
>37.4 - 38.5 °C (>99.4 - 101.3°F)	Elevated
>38.5 - 42.2 °C (>101.3 - 108.0°F)	Possible high fever

*Fever above 40°C (104°F) require immediate medical attention especially in young children.

Stethoscope setup

- a. Twist Cap off
- b. Push down on child proof lock and twist open.
- c. Insert battery
- d. Re-assemble


Taking Recordings – Stethoscope

Step 1) Make sure stethoscope is connected to the smartphone, the LED indicator will turn on briefly.

Note: If LED on the stethoscope does not light up when connected to the smartphone, there may be a problem with your 3.5mm headphone socket on your phone. Commonly, this is caused by debris accumulated in the jack. Please contact your phone manufacturer if this problem persists.

Step 2) Open heart or lung recordings by following on-screen instructions in the app. The LED indicator will turn on.

Step 3) Position the stethoscope against locations specified by in-app instructions. The diaphragm and blue o-ring should be pressed against contact surface. For more information, see in-app tutorials.

Hold stethoscope with fingers around the grooves on the side. Keep hands steady to achieve best results.

Step 4) Initiate recording in the app. Apply enough pressure to keep stethoscope from moving position, but not too much pressure as to cause discomfort. During recording in app, keep noise and movements to a minimum.

Note: Recordings against the skin will provide the best sound quality recording. If you do record through clothing, avoid doing so through thick garments.

Step 5) Playback the recording in the app, through headphones.

Note: The CliniCloud Stethoscope calculates heart rate by counting the number of beats in a heart sound recording. This measurement may be affected by background noise, or movements of the stethoscope during recording. For the most accurate heart rate or monitoring of clinical vitals, we recommend the use of an electrocardiogram (ECG).

Step 6) To turn off device, simply unplug the stethoscope from the phone.

Note: Interpreting stethoscope recordings require professional training and self-diagnoses is not recommended. CliniCloud lets you share these recordings via email or through our integration with telemedicine providers, such as Doctor On Demand.

Sending Recordings

Stethoscope recordings of the heart and lung may contain information that can be used to diagnose and monitor severity of diseases, like asthma or bronchiolitis. Please check in-app instructions for how to send recordings to your physician or another party. Note that some services (including on demand doctor review of recordings) may not be available in your region. Please check <http://clinicloud.com> for the latest information.

Device maintenance

In order to extend the lifetime of the devices and components, please observe the following conditions:

- Do not expose the devices to high temperatures, direct sunlight, high humidity, highly dusted areas or strong vibrations.
- Do not disassemble, repair or modify the devices. If you have any questions, please contact the manufacturer's customer service.
- Avoid dropping, slamming or throwing devices
- Keep devices dry at all times
- Batteries should be removed from devices, during extended periods of disuse

Cleaning instructions

To clean the entire device, please use a soft cloth slightly moistened with a mild soapy solution. Under no circumstances may liquid enter the device. Do not use the device again until it is completely dry.

Useful life

The useful life of the devices is 5 years.

Support

Warranty : 1 Year Limited Warranty

Manufacturer: StethoCloud Pty. Ltd.

(Subsidiary of CliniCloud Inc.)

Check <http://clinicloud.com> for further details

Email : support@clinicloud.com

Address : Level 5, 141 Flinders Lane, Melbourne 3000, Australia

**European Authorized
Representative :**

Emergo Europe
Molenstraat 15
2513 BH The Hague
The Netherlands

© 2015 CliniCloud, Inc. All Rights Reserved. CliniCloud is a pending trademark of CliniCloud, Inc. CliniCloud products are protected by US and foreign patents.

Visit <http://clinicloud.com> for more information.

Standards

Declaration of Conformity

We CliniCloud, Inc.*

Address:

350 Townsend #758

San Francisco, 94107

California, USA

+1 (415) 801 3283

Declare under our sole responsibility that the product, The CliniCloud Thermometer, complies with part 15 of the FCC rules. Operation is subject to the following two conditions: **(1)** this device may not cause harmful interference, and **(2)** this device must accept any interference received, including interference that may cause undesired operation.

FCC Caution: changes / modifications not approved by the responsible party could void the user's authority to operate the devices.

FCC ID: 2AE7C-SPL1024

This device (CliniCloud Thermometer) complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions: **(1)** this device may not cause interference, and **(2)** this device must accept any interference, including interference that may cause undesired operation of the device.

IC : 20418-SPL1024

This CliniCloud infrared thermometer also meets requirements established in ASTM Standard (E1965-98). Full responsibility for the conformance of this product to the standard is assumed by CliniCloud Inc.

Please kindly be aware that it is the consumer's responsibility to properly dispose and recycle devices and accompanying components. The devices are considered electronic waste and should be disposed of at your local electronic equipment collection facility. Please dispose of used batteries according to your local environmental protection requirements.

*CliniCloud Inc is the parent company of the manufacturer, StethoCloud Pty. Ltd.

Troubleshooting

Stethoscope

Problem	Solution
My stethoscope isn't working.	<ol style="list-style-type: none">1. Make sure the battery is installed properly.2. Make sure the stethoscope is connected to the phone. The LED should light up briefly once connected. Note: If LED doesn't light up, there's likely an issue with your phone's headphone port. Contact your phone manufacturer for more details.3. Make sure you have downloaded and logged into the app.4. Once stethoscope is plugged in, try turning up the headphone volume.5. Try re-plugging the stethoscope.
I can't capture my heart or lung sounds.	Ensure your stethoscope is working and connected, and try again in a quiet environment. Please place stethoscope as instructed in the app user interface.
The LED light is flashing more quickly than usual.	This is a low battery signal, please change batteries at your earliest convenience.
I can't use other apps with the CliniCloud Stethoscope.	The device is only designed to be used with the CliniCloud app.

Thermometer

Problem	Solution
My thermometer isn't working.	<ol style="list-style-type: none">1. Make sure the battery is installed properly.2. Make sure Bluetooth is turned on.3. Make sure you have downloaded and logged into the app.4. Enter the Temperature applet – the thermometer LED should light up.5. During the scan, the thermometer LED should flash intermittently.6. If this doesn't work, try resetting the thermometer by reinstalling the batteries.
I'm getting an inaccurate temperature.	Ensure you are scanning the forehead correctly. The thermometer takes multiple readings during the scan, and displays only the highest reading.
I have multiple thermometers.	Great! If they are both within your Bluetooth range, you will be prompted in the app which thermometer you wish to pair with. Follow on-screen UI to choose the right one.
I can't pair with the thermometer.	You don't need to pair to use the thermometer. Do not attempt to pair with the device.
I can't use other apps with the CliniCloud Thermometer.	The device is only designed to be used with the CliniCloud app.

Important Warnings

- Prolonged high fever requires medical attention especially in young Children.
- Wheezing and respiratory issues can progress quickly and be fatal the CliniCloud stethoscope is not made for self diagnosis and results should be referred to a professional. 911 (or local emergency number) should be called in the event of any respiratory event that causes concern.
- The thermometer needs to be in the room in which the measurement is taken for at least 15 minutes before use. Do not place thermometer on heated surface between uses, as the thermometer needs to be adjusted to room temperature for accurate readings.
- The CliniCloud thermometer is only designed for measuring points on the human body stated in the instructions for use. If measuring temperature on a child in a radiant warmer, use oral or rectal thermometers instead.
- The devices are only intended for the purpose stated in these instructions for use. Children must not be allowed to use the devices. Medical products are not toys.
- The devices have been designed for practical use, but is not a substitute for a visit to the doctor. Should you have any questions about using the device, please contact either your retailer or Customer Services.
- Packaging materials are a deadly hazard for children and can cause suffocation. Remove all packaging materials immediately and keep them away from Children at all times.
- Swallowing batteries and/or battery fluid can be extremely dangerous. Keep the battery out of reach of Children. Should a child or any other person swallow a battery and/or battery fluid, please call 911 or your local emergency number immediately.
- The battery cap on the stethoscope was designed with a child safety lock. Please do not remove coin cell battery from the stethoscope between uses to avoid children from accidentally swallowing battery.
- Should battery fluid leak and come into contact with your eyes or skin, immediately rinse with plenty of clean water. Please call 911 or your local emergency number immediately.
- Do not take recordings if devices are exposed to strong electromagnetic field or static electricity.
- Thermometer device performance may be adversely affected if:
 - Operated outside of the specified
 - Operated outside of the specified operating temperature and humidity ranges.
 - Stored outside of the manufacturer specified ambient temperature and humidity ranges.
 - Subject to mechanical shock

Limited Warranty

This Warranty for the CliniCloud thermometer and stethoscope (the "Devices") and is provided by CliniCloud Inc,
350 Townsend #758
San Francisco, 94107
California, USA

The provider of this warranty is sometimes referred to herein as "we". We warrant the Devices against defects in materials and workmanship under ordinary consumer use for one year from the date of original retail purchase. During this warranty period, if a defect arises in the Devices, and you follow the instructions for returning the Devices, we will, at our option and to the extent permitted by law, either (i) repair the Devices using either new or refurbished parts, (ii) replace the Devices with a new or refurbished Devices, or (iii) refund to you all or part of the purchase price of the Devices. This limited warranty applies to any repair, replacement part or replacement Devices for the remainder of the original warranty period or for ninety days, whichever period is longer. All replaced parts and Devices for which a refund is given shall become our property. This limited warranty applies only to hardware components of the Devices that are not subject to accident, misuse, neglect, fire or other external causes, alterations, repair, or commercial use.

Instructions

For specific instructions about how to obtain warranty service for your Device, please contact Customer Service using the contact information provided in this guide. In general, you will need to deliver your Devices in either its original packaging or in equally protective packaging to the address specified by Customer Service. Before you deliver your Devices for warranty service, it is your responsibility to back up any data, software, or other materials you may have stored or preserved on your Devices. It is possible that such data, software or other materials will be lost or reformatted during service, and we will not be responsible for any such damage or loss.

Limitations

TO THE EXTENT PERMITTED BY LAW, THE WARRANTY AND REMEDIES SET FORTH ABOVE ARE EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES AND REMEDIES, AND WE SPECIFICALLY DISCLAIM ALL STATUTORY OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND AGAINST HIDDEN OR LATENT DEFECTS. IF WE CANNOT LAWFULLY DISCLAIM STATUTORY OR IMPLIED WARRANTIES, THEN TO THE EXTENT PERMITTED BY LAW, ALL SUCH WARRANTIES SHALL BE LIMITED IN DURATION TO THE DURATION OF THIS EXPRESS LIMITED WARRANTY AND TO REPAIR OR REPLACEMENT SERVICE. SOME JURISDICTIONS DO NOT ALLOW LIMITATIONS ON HOW LONG A STATUTORY OR IMPLIED WARRANTY LASTS, SO THE ABOVE LIMITATION MAY NOT APPLY TO YOU. WE ARE NOT RESPONSIBLE FOR DIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM ANY BREACH OF WARRANTY OR UNDER ANY OTHER LEGAL THEORY. IN SOME JURISDICTIONS THE FOREGOING LIMITATION DOES NOT APPLY TO DEATH OR PERSONAL INJURY CLAIMS, OR ANY STATUTORY LIABILITY FOR INTENTIONAL AND GROSS NEGLIGENT ACTS AND/OR OMISSIONS, SO THE ABOVE EXCLUSION OR LIMITATION MAY NOT APPLY TO YOU. SOME JURISDICTIONS DO NOT ALLOW THE EXCLUSION OR LIMITATION OF DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE ABOVE EXCLUSION OR LIMITATION MAY NOT APPLY TO YOU.

This limited warranty gives you specific legal rights. You may have additional rights which vary from jurisdiction to jurisdiction, and this limited warranty does not affect those rights.

Technical information

Thermometer

Measurement method	Non-contact infrared measurement
Measurement units	Fahrenheit (°F) and Celsius (°C)
Recommended operating conditions	59°F to 104°F (16°C to 40°C) with a relative humidity of ≤85%
Storage conditions	-4°F to 122°F (-20°C to 50°C) with a relative humidity of ≤85%, 700 hPa to 1060 hPa
Measurement distance	0" to 2" (0 to 5cm) from the measuring point
Measurement range	90.0°F to 109.9°F (32.2°C to 43.3°C)
Measurement accuracy	± 0.4°F (±0.2°C)
Energy saving functions	Device automatically returns to low power advertising mode after 2 minutes
Dimensions (Width x depth x height)	approx. 0.98" x 1.42" x 4.53" (approx. 25mm x 36 mm x 155mm)
Connectivity	Bluetooth Low Energy
Power supply	3V, 2 x AAA
Expected battery life	1 year, with single daily use

Smartphone minimum requirement:

iOS

Hardware	Any iOS device with Bluetooth 4.0 compatibility, i.e. iPhone 4S and newer, iPod Touch 5th Gen and newer.
Operating system	iOS 8.0 and higher

Android

Hardware

Any Android smartphone with Min 512mb of storage space, and Bluetooth 4.0 compatibility.

Operating System

Android 4.3 (Jelly Bean) and higher

Stethoscope

Basic functions

Heart sound recording, Lung sound recording

Operating conditions

50°F to 104°F (10°C to 40°C) with a relative humidity of ≤85%

Storage conditions

-4°F to 122°F (-20°C to 50°C) with a relative humidity of ≤85%, 700 hPa to 1060 hPa

Recording quality

44.1 kHz 16 bit Linear PCM via CliniCloud App (when recorded on iPhone 6)

Energy saving functions

Device goes to sleep when not connected to smartphone

Dimensions (Diameter x depth)

approx. 1.42" x 0.98" (approx. 36 mm x 25 mm)

Connectivity

micro-USB to 3.5mm jack cord

Power supply

3V Coin cell (CR2032)

Expected battery life

1 year+ on standby, 10+ hours of recording time

Smartphone minimum requirement:

iOS

Hardware

Any iOS device with Bluetooth 4.0 compatibility, i.e. iPhone 4S and newer, iPod Touch 5th Gen and newer.

Operating system

iOS 8.0 and higher

Android

Hardware	Any Android smartphone with Min 512mb of storage space, and Bluetooth 4.0 compatibility.
Operating System	Android 4.3 (Jelly Bean) and higher

Symbols Used

	European conformity mark
	WEEE – Properly dispose of Electronic Waste
	Consult Instructions for Use / Manual
	Manufacturer
	Temperature Limits (Storage)
	Batch code / Lot number
	Date of manufacture
	Keep dry
	Authorized European Representative
	Type BF applied part
	Catalogue number

Stethoscope:

Guidance and manufacturer’s declaration – electromagnetic emission – for all EQUIPMENT and SYSTEMS

Guidance and manufacturer’s declaration – electromagnetic emissions

The model SPL999 is intended for use in the electromagnetic environment specified below. The customer or the user.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Model SPL999 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions	Class B	The Model SPL999 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations / flicker	N/A	

Guidance and manufacturer’s declaration – electromagnetic emission – for all EQUIPMENT and SYSTEMS

Guidance and manufacturer’s declaration – electromagnetic emissions

The model SPL999 is intended for use in the electromagnetic environment specified below. The customer or the user.

Immunity test	EN 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	N/A	
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	N/A	
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	N/A	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE— U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer’s declaration – electromagnetic immunity – for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer’s declaration—electromagnetic immunity

The model SPL999 is intended for use in the electromagnetic environment specified below. The customer or the user of the model SPL999 should assure that it is used in such an environment.

Immunity test	EN 60601 test	Compliance level	Electromagnetic environment – guidance
---------------	---------------	------------------	--

Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A	<p>Portable and mobile RF communications equipment should be used no closer to any part of the SPL999, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ <p>$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz</p> $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1—At 80 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2—These guidelines may not apply in all situations.
 Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people

- a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model SPL999 is used exceeds the applicable RF compliance level above, the model SPL999 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the model SPL999.
- b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM – for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the model SPL999

The model SPL999 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the model SPL999 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model SPL999 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{\sqrt{f}} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{\sqrt{f}} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{\sqrt{f}} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1—At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2—These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Description of wireless technology

The CliniCloud Thermometer uses Bluetooth LE technology to transmit temperature from the thermometer to the smartphone app. Please note the following information:

- To guarantee safe and effective operation of the device, the thermometer should be located within 15 meters (line-of-sight) and ideally in the immediate vicinity of the smartphone for the entire duration of the recording.
- Recordings made using the CliniCloud non-contact thermometer are transferred directly to the smartphone app and the thermometer itself does not store user data. The device implements the standard Bluetooth GATT based health thermometer profile. Wireless security measures are included as part of the Bluetooth standard.
- If problems occur and the CliniCloud non-contact thermometer fails to establish/maintain connection with the smart phone, try to decrease distance between the device and the smartphone, and move away from strong EMI sources such as wireless routers. If problems persist, please contact support@clinicloud.com.

Thermometer:

Guidance and manufacturer’s declaration – electromagnetic emission – for all EQUIPMENT and SYSTEMS

Guidance and manufacturer’s declaration – electromagnetic emissions

The model SPL1024 is intended for use in the electromagnetic environment specified below. The customer or the user

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Model SPL1024 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Model SPL1024 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	N/A	

Guidance and manufacturer’s declaration – electromagnetic emission – for all EQUIPMENT and SYSTEMS

Guidance and manufacturer’s declaration—electromagnetic immunity

The model SPL1024 is intended for use in the electromagnetic environment specified below. The customer or the user of the model SPL1024 should assure that it is used in such an environment.

Immunity test	EN 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	N/A	
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	N/A	
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	N/A	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE— U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer’s declaration – electromagnetic immunity – for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer’s declaration—electromagnetic immunity

The model SPL1024 is intended for use in the electromagnetic environment specified below. The customer or the user of the model SPL1024 should assure that it is used in such an environment.

Immunity test	EN 60601 test level	Compliance level	Electromagnetic environment – guidance
---------------	---------------------	------------------	--

Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the model SPL1024, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance
			$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \text{ 800 MHz to 2.5 GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol:</p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1—At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2—These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people

- a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model SPL1024 is used exceeds the applicable RF compliance level above, the model SPL1024 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the model SPL1024.
- b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM – for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the model SPL1024

The model SPL1024 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the model SPL1024 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model SPL1024 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{f} \right] P$	80 MHz to 800 MHz $d = \left[\frac{3.5}{f} \right] P$	800 MHz to 2.5 GHz $d = \left[\frac{7}{f} \right] P$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	100
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1—At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2—These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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